

QUARTERLY COMPOSITING OF AIRNET SAMPLES

Purpose This Meteorology and Air Quality Group (MAQ) procedure describes the preparation of composites of biweekly AIRNET samples for the purpose of determining alpha-emitting nuclides and selected inorganics on a quarterly basis.

Scope This procedure applies to the preparation and submittal of composited samples of filter papers used to collect airborne particulates as part of the AIRNET monitoring program.

In this procedure This procedure addresses the following major topics:

Topic	See Page
General Information About This Procedure	2
Who Requires Training to This Procedure?	2
Preparation for Quarterly Sample Compositing	4
Compositing Samples at the Laboratory	6
Records Resulting from This Procedure	9

Hazard Control Plan The hazard evaluation associated with this work is documented in Attachment 1: Initial risk = **minimal**. Residual risk = **minimal**. Work permits required: none.
First authorization review date is one year from group leader signature below; subsequent authorizations are on file in group office.

Signatures

Prepared by: Ernie Gladney, MAQ	Date: <u>6/17/2002</u>
Approved by: Craig Eberhart, Air Quality Monitoring Project Leader	Date: <u>6/17/2002</u>
Approved by: Terry Morgan, QA Officer	Date: <u>6/17/02</u>
Work authorized by: Jean Dewart, MAQ Acting Group Leader	Date: <u>6/18/02</u>

07/11/02

CONTROLLED DOCUMENT

This copy is uncontrolled if no signatures are present or if the copy number stamp is black. Users are responsible for ensuring they work to the latest approved revision.

General information about this procedure

Attachments This procedure has the following attachments:

Number	Attachment Title	No. of pages
1	Hazard Control Plan	2

History of revision

This table lists the revision history and effective dates of this procedure.

Revision	Date	Description Of Changes
0	3/13/2001	New document created from a chapter of ESH-17-202, R12.
1	7/2/02	Updated details about station 90 handling, cutting filters, and changes to requested analytes.

Who requires training to this procedure?

The following personnel require training before implementing this procedure:

- MAQ personnel assigned to composite AIRNET samples

Annual retraining is required and will be by self-study (“reading”) training.

Training method

The training method for this procedure is **on-the-job** training by a previously trained individual and is documented in accordance with the procedure for training (MAQ-024).

Prerequisites

In addition to training to this procedure, the following training is also required prior to performing this procedure:

- Appropriate Radiation Worker training (as required by the location where procedure is performed)

General information, continued

**Definitions
specific to this
procedure**

None.

References

The following documents are referenced in this procedure:

- MAQ-024, “Personnel Training”
 - MAQ-202, “Environmental Sampling of Airborne Particulate Radionuclides”
 - MAQ-208, “Evaluation of Biweekly AIRNET Data and Calculated Air Concentrations”
 - MAQ-AIRNET, “Sampling and Analysis Plan for the Radiological Air Sampling Network (AIRNET)”
-

Note

Actions specified within this procedure, unless preceded with “should” or “may,” are to be considered mandatory guidance (i.e., “shall”).

Preparation for quarterly sample compositing

Overview

AIRNET filter samples are collected every two weeks according to procedure MAQ-202. The uncut biweekly samples are shipped to the analytical laboratory (currently in Grand Junction) for instrumental analysis. Every quarter, the filters from the past calendar quarter are composited into one sample per each AIRNET site for analysis for additional radionuclides and elements. The compositing is done at the laboratory in conjunction with one of that laboratory's staff. Thus, two people commonly perform this process: one to prepare the paperwork and oversee the work, and the second person to handle the samples. See MAQ-AIRNET for more details.

Preparing for quarterly composite

Compositing of the quarterly samples is performed at the analytical laboratory at the end of the quarter after the final set of instrumental analyses is completed. Perform the following steps to prepare needed paperwork for compositing.

Step	Action
1	In the MS Access AIRNET database AIRNET_PeriodID table, there is a quarterly composite field, "QuartComp", that will indicate which biweekly period IDs will be included in each specific composite. This field is populated during each calendar year's start-up activities. NOTE: Samples deployed before Christmas and collected in early January of the following year will be included in the first quarterly composite of the new year.
2	Ensure the field data V&V for filters is complete and the technical review has been performed according to MAQ-208. If the memo documenting the technical review is not finalized, obtain e-mail or other written approval from the Air Quality Monitoring Project Leader. NOTE: We elect to take the risk that a few samples may be reclassified upon Health Physics review, and will NOT use the gross beta as an indicator for compositing decisions and may composite prior to completion of that review.

Steps continued on next page.

Preparation for quarterly sample compositing, continued

Step	Action
3	Using the MS Access AIRNET database (see the AIRNET Database User's Guide), produce the composite reports to ensure that no sample that has been rejected for field data reasons gets inadvertently included in a composite. This is accomplished by starting on the AIRNET Database Form "MAIN Switchboard" and selecting the "Quarterly Composite" button. The subsequent forms will walk you through the process of producing a complete checklist for assembling the composite, printing labels, editing the cover letter, producing a chain-of-custody document for returning samples, etc. Be sure to take a copy of all documents produced via this database form to the compositing site.
4	Prepare pre-printed labels for the composited samples using the button for that purpose on the AIRNET Database "Quarterly Composite" Form. Be sure to take these with you to the compositing site.
5	Prepare a letter to the analytical laboratory by editing the existing template (via the database form) that requests the analyses on the quarterly composite samples. The chemistry data coordinator maintains current analysis information and it is also documented in the Sites_MasterLocation table within the AIRNET database.

Compositing samples at the laboratory

Compositing the samples

After the filters have been collected, shipped, and analyzed for one quarter, they are ready to be “composited” and submitted for the analyses of various radionuclides by radiochemical alpha spectroscopy. This process will be performed at the analytical laboratory under the direction of MAQ personnel.

Supplies needed

Collect the materials and tools listed below.

- gloves (powderless)
- scissors
- metal sample cans
- Kimwipes
- 2 sets of labels (e.g., “00Q1.nn”)
- tweezers
- Copy of MS Access reports or forms for “Quarterly Composites”
- bench paper
- permanent marker
- scotch tape
- cleaning solution

Steps to composite the samples for analysis

To composite the samples, perform the following steps:

Step	Action
1	Ensure all required training is complete and current for working in the analytical laboratory. A Grand Junction Analytical Laboratory (GJAL) person is required to perform the actual filter cutting operation with participating MAQ staff to help track the filters, record the samples composited on the Access report printout and take legal responsibility for the overall compositing process.
2	The GJAL staff will collect all filters (in individual metal cans) for the calendar quarter to be submitted. Refer to the Compositing Checklist generated in step three on previous page.
3	For each sample station, label (or re-label) one metal sample can (or other appropriate containers) with the station I.D., in the format yyQn.ss, where yy is the year, n is the quarter number, and ss is the station number, using the preprinted labels generated from the AIRNET database. These labels will be generated from the AIRNET database and supplied by the participating MAQ staff member.
4	Put all sample cans (each containing an individual filter) from each station in an individual stack.

Steps continued on next page.

Compositing samples at the laboratory, continued

Step	Action
5	Refer to the Compositing Checklist generated above and remove all samples marked with an “R”. Rejected samples will be analyzed individually. For these rejected samples, label a new can (or other container), put an ‘X’ after the sample ID on the sample’s envelope, and enter the sample ID on the compositing form (e.g., 01Q1.14X). A preprinted label should be available for each rejected sample.
6	Refer to the Compositing Checklist and, with the assistance of the second person, ensure the stack contains all the period IDs listed on the checklist. For each sample that goes in a composite, mark on the Compositing Checklist that the samples are present and were included in this composite.
7	Cover the bench with bench paper and place a Kimwipe over the immediate work area. Put on the gloves. Select one stack (representing a single station).
8	Place a labeled sample composite can open, face-up in the covered work area. Remove samples one at a time from their original counting containers, and holding each over the composite can with the tweezers, cut it in half, letting the loose half drop gently into the composite can. Place the remaining half (held in the tweezers) gently back into its original counting can. Note that as a result of changes to the biweekly analysis requirements, that some of the stations will have only one-half filter remaining for use in the composite. In these cases, simply combine all the remaining half filters WITHOUT any further cutting.
9	If there is remaining loose particulate material in any of the vacated individual sample cans, manually apportion it uniformly between the new composite cans.
10	Discard any vacated individual sample cans. Repeat steps 8 and 9 for the remaining filters in this quarter for this site. Close the completed composite sample can and set aside for resubmission to the Sample Management organization at the analytical laboratory. Since the filters have never left their custody, all chain-of-custody documentation is handled internal to the analytical laboratory.
11	Put a new Kimwipe on the work area. Wash both scissors and tweezers. Select the next stack (representing a single station). Repeat steps 8 through 10.
12	After the last composite is prepared, both the participating MAQ staff and the analytical laboratory staff must sign and date the Compositing Checklist. Make a copy, and leave one copy at the analytical laboratory. Leave the original of the cover letter requesting the specific analyses to be performed. Ensure the original is returned and becomes part of the records retained at MAQ.

Steps continued on next page.

Compositing samples at the laboratory, continued

Step	Action
13	Package all remaining sample half-filters for shipment back to Los Alamos. Initiate the C-of-C form that was prepared from the database and receive custody of the sample remains from the analytical laboratory for return to MAQ. Upon arrival in Los Alamos, relinquish custody to the field team at our TA-54 site for sample archiving.

Records resulting from this procedure

Records

The following records generated as a result of this procedure are to be submitted **within 1 week** as records to the records coordinator:

- Completed and signed Compositing Checklist. This will generally be a part of a larger internal memo that documents how the quaterly air volumes were calculated.

HAZARD CONTROL PLAN

1. The work to be performed is described in this procedure.

“Quarterly Compositing of AIRNET Samples”

2. Describe potential hazards associated with the work (use continuation page if needed).

Contamination from handling filters with potential particulate activity.

3. For each hazard, list the likelihood and severity, and the resulting initial risk level (before any work controls are applied, as determined according to LIR300-00-01.0, section 7.2)

Improbable / Negligible = Minimal

Overall *initial* risk: ☒ Minimal ☐ Low ☐ Medium ☐ High

4. Applicable Laboratory, facility, or activity operational requirements directly related to the work:

☐ None ☒ List:

Work Permits required? ☒ No ☐ List:

LIR-402-706-01 “Personnel Dosimetry”

HAZARD CONTROL PLAN, continued

5. Describe how the hazards listed above will be mitigated (e.g., safety equipment, administrative controls, etc.):

Use of appropriate PPE (gloves, lab coat), work done in a controlled area. Appropriate radiological monitor available, RCT phone number posted by phone in case of contamination.

6. Knowledge, skills, abilities, and training necessary to safely perform this work (check one or both):



Group-level orientation (per MAQ-032) and training to this procedure.



Other → See training prerequisites on procedure page 3. Any additional describe here:

7. Any wastes and/or residual materials? (check one) ☒ None ☐ List:
(Used filter media or analytical residue is disposed by the analytical laboratory.)

8. Considering the administrative and engineering controls to be used, the *residual* risk level (as determined according to LIR300-00-01.0, section 7.3.3) is (check one):



Minimal



Low



Medium (requires approval by Division Director)

9. Emergency actions to take in event of control failures or abnormal operation (check one):



None



List:

A: Notify ESH-1 RCT for area or other appropriate rad control personnel.

B: Use decontamination material such as “Fantastic” cleaner, paper towel, gloves.

C: Notify MAQ supervisor.

Signature of preparer of this HCP: This HCP was prepared by a knowledgeable individual and reviewed in accordance with requirements in LIR 300-00-01 and LIR 300-00-02.

Preparer(s) signature(s)

Name(s) (print)

/Position

Date

Signature by group leader on procedure title page signifies authorization to perform work for personnel properly trained to this procedure. This authorization will be renewed annually and documented in MAQ records.

Controlled copies are considered authorized. Work will be performed to controlled copies only. This plan and procedure will be revised according to MAQ-022 and distributed according to MAQ-030.

